

These are approximate flow rate / FiO₂

NASAL CANNULA



O ₂ flow (l/min)		FiO ₂ (%)
1	→	24
2	→	28
3	→	32
4	→	36
5	→	40

FACE MASKS



VENTURI MASK

VALVE	O ₂ flow (l/min)		FiO ₂ (%)
BLUE	2	→	24
WHITE	4	→	28
ORANGE	6	→	31
YELLOW	8	→	35
RED	10	→	40
GREEN	15	→	60



SIMPLE FACE MASK

O ₂ flow (l/min)		FiO ₂ (%)
5	→	40
6-7	→	50
7-8	→	60



NON-REBREATHER RESERVOIR

Ensure reservoir is filled with oxygen prior to use

O ₂ flow (l/min)		FiO ₂ (%)
15	→	85-90